


## SUGGESTED SET LIST

- SECOND DEGREE TEAR - Injury to the perineum involving perineal muscles but not involving the anal sphincter.
- THIRD DEGREE TEAR - Injury to the perineum involving the anal sphincter complex.
- FOURTH DEGREE TEAR - Injury to the perineum involving the anal sphincter complex and rectal mucosa.



NB. We do not supply consumables (such as gauze, plastic holloware and solutions)

This brochure was produced with the kind assistance of Mr Gelman, (Consultant Urogynaecologist) and Lead Practitioner Carolyn Elson at St. Peter's Hospital in Chertsey, England.

For samples or further information on any product, please contact our sales team on (020) 7511 9960.

**Timesco**<sup>TM</sup>  
SURGICAL + MEDICAL  
INSTRUMENTS FOR LIFE<sup>TM</sup>

Timesco of London Ltd, 1 Knights Road, London E16 2AT England  
Telephone: +44 (020) 7511 9960 Fax: +44 (020) 7511 7888 Email: [surgical@timesco.com](mailto:surgical@timesco.com)  
[www.timesco.com](http://www.timesco.com)

E & O.E. Copyright April 2004. Ref C.S. 005/1  
'Timesco' and 'Instruments for Life' are registered trademarks of Timesco of London. All rights reserved.




Please visit our website  
at [www.timesco.com](http://www.timesco.com)



THIRD AND FOURTH DEGREE  
OBSTETRIC TEAR SET INSTRUMENTS

**Timesco**<sup>TM</sup>  
SURGICAL + MEDICAL  
INSTRUMENTS FOR LIFE<sup>TM</sup>

# THIRD AND FOURTH DEGREE TEAR SET INSTRUMENTS



*Our Officers Forceps have small 3x4 teeth which provide secure grip with minimal trauma.*


Our Obstetric and Gynaecology Instruments are vital to our extensive range of theatre surgical Instruments. With the kind support and guidance of the Urogynae team at St. Peter's Hospital in Chertsey, we have produced this brochure as a guide for instruments used in this vital procedure. Each instrument has been carefully selected to compile the set listed on the back page.

Although we recommend these instruments for this procedure, other instruments may be supplemented from our extensive range, based on each Surgeons preference.

All our reusable theatre instruments are produced from the finest surgical stainless steels to exacting tolerances to deliver precise clinical performance. All instruments are fully autoclavable and covered by our unrivalled lifetime guarantee.




14.070.10 CRILE ARTERY FORCEPS 5.5" (14CM)




16.600.20 ALLIS TISSUE FORCEPS 3X4 TEETH 6" (15.25CM)



12.500.10 MAYO SCISSORS STRAIGHT 6.5" (16.5CM)



12.965.20 MCINDOE SCISSORS CURVED 7" (17.75CM) BLUNT/BLUNT T.C.



16.100.40 RAMPLEY SPONGE HOLDING FORCEPS 9.5" (24.25CM)



17.060.10 CRILE MURRAY NEEDLE HOLDER 6" (15.25CM)



11.550.10 BP SCALPEL HOLDER No.3 4.5" (11.5CM)



15.675.15 WEISLANDER SELF RETAINING RETRACTOR 3X4 PRONGS BLUNT 5.5" (14CM)



13.325.10 OFFICER DISSECTING FORCEPS 4X5 TEETH SERRATED 6" (15.25CM)



15.425.10 LANDON VAGINAL RETRACTOR 8" (20.25CM)  
1"X3.5" (2.54X9CM)



57.575.30  
SIMS VAGINAL SPECULUM LARGE  
DOUBLE ENDED



57.575.40  
SIMS VAGINAL SPECULUM MEDIUM  
DOUBLE ENDED



*We recommend the use of a blunt jawed Weislander Retractor to minimise tissue trauma to the labia.*